

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 124 255

**Pellet-burettes, 25 ml, amber glass, AS, conf.cert.**




**Batch:** 43-20

**Serial:** -

**Measuring Nr.:** 59623

**Batch size:** 32

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	24,9795 ml	14,9890 ml	4,9893 ml
Deviation:	0,0021 ml	0,0035 ml	0,0035 ml
Precision (CV):	0,0083 %	0,0230 %	0,0702 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

4

**Date:**

06.10.2022

**Balance:**

Sartorius -  
MSE224S-0CE  
-DU

