

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 124 505

**Pellet-burettes, 50 ml, amber glass, AS, conf. cert.**




**Batch:** 25-20

**Serial:**

**Measuring Nr.:** 47060

**Batch size:** 10

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	49,9918 ml	29,9836 ml	10,0107 ml
Deviation:	0,0175 ml	0,0075 ml	0,0108 ml
Precision (CV):	0,0351 %	0,0249 %	0,1082 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:** 77

**Date:** 28.05.2021

**Balance:** Sartorius -  
LE244S-OCE

DE-M