

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 124 508

**Pellet-burettes, 50 ml, amber glass, AS, conf. cert.**




**Batch:** 18-21

**Serial:**

**Measuring Nr.:** 49639

**Batch size:** 15

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0012 ml	29,9618 ml	10,0089 ml
Deviation:	0,0072 ml	0,0013 ml	0,0086 ml
Precision (CV):	0,0144 %	0,0044 %	0,0864 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	07.10.2021
Balance:	Sartorius - MSE224S-0CE -DU

