

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 124 508

Pellet-burettes, 50 ml, amber glass, AS, conf. cert.




Batch: 43-21

Serial:

Measuring Nr. 56231

Batch size: 10

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0015 ml	29,9852 ml	10,0285 ml
Deviation:	0,0125 ml	0,0167 ml	0,0030 ml
Precision (CV):	0,0250 %	0,0555 %	0,0302 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

25.05.2022

Balance:

Sartorius -
LE244S-0CE

DE-M