

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 003 2 BL

Grad. cylinder, 25:0,2 ml, blue grad., A, conf. cert.




Batch: 30-24

Serial:

Measuring Nr. 68427

Batch size: 20

Volume:	5,0000 ml	15,0000 ml	25,0000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	5,0395 ml	15,0274 ml	25,0254 ml
Deviation:	0,0386 ml	0,0275 ml	0,0457 ml
Precision (CV):	0,7714 %	0,1835 %	0,1830 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	06.08.2024
Balance:	Sartorius - MSE224S-0CE -DU

