

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 005 BL

**Graduated cylinders, 50:1 ml, blue grad., A, conf. cert.**




**Batch:** 51-24

**Serial:**

**Measuring Nr.:** 71241

**Batch size:** 400

Volume:	10,000 ml	30,000 ml	50,000 ml
<b>Tolerance:</b>	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
<b>Average value:</b>	9,9021 ml	30,0234 ml	50,0791 ml
<b>Deviation:</b>	0,0874 ml	0,0675 ml	0,0695 ml
<b>Precision (CV):</b>	0,8741 %	0,2250 %	0,1389 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	18.12.2024
<b>Balance:</b>	Sartorius - MSE224S-0CE -DU

