

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 010 BL

**Grad. cylinder, 100:1,0 ml, blue grad., A, conf. cert.**




**Batch:** 39-24

**Serial:**

**Measuring Nr.** 69416

**Batch size:** 500

Volume:	20,000 ml	60,000 ml	100,000 ml
<b>Tolerance:</b>	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
<b>Average value:</b>	20,0530 ml	60,0276 ml	99,9902 ml
<b>Deviation:</b>	0,1167 ml	0,1060 ml	0,0982 ml
<b>Precision (CV):</b>	0,5834 %	0,1766 %	0,0982 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	07.10.2024
<b>Balance:</b>	Sartorius - LA 1200 S

