

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 010

Graduated cylinder, 100:1,0 ml, A, conf. cert.




Batch: 39-24

Serial:

Measuring Nr.: 69452

Batch size: 1200

Volume:	20,000 ml	60,000 ml	100,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	20,1296 ml	60,0186 ml	100,0194 ml
Deviation:	0,1062 ml	0,0799 ml	0,0905 ml
Precision (CV):	0,5310 %	0,1331 %	0,0905 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	11.10.2024
Balance:	Sartorius - MSE224S-0CE -DU

