

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 050 USP

Grad. cylinder, 500:5 ml, USP, A, conf. cert.




Batch: 27-24

Serial:

Measuring Nr.: 67845

Batch size: 280

Volume:	100,000 ml	300,000 ml	500,000 ml
Tolerance:	± 2,0000 ml	± 2,0000 ml	± 2,0000 ml
Average value:	99,6681 ml	300,1620 ml	500,2336 ml
Deviation:	0,3581 ml	0,3328 ml	0,3123 ml
Precision (CV):	0,3581 %	0,1109 %	0,0625 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	04.07.2024
Balance:	Sartorius - LA 1200 S

