

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 100 BL

**Grad. cylinder, 1000:10 ml, blue grad., A, conf. cert.**




**Batch:** 31-24

**Serial:**

**Measuring Nr.** 68530

**Batch size:** 280

Volume:	200,000 ml	600,000 ml	1000,000 ml
<b>Tolerance:</b>	± 5,0000 ml	± 5,0000 ml	± 5,0000 ml
<b>Average value:</b>	200,3051 ml	599,8426 ml	999,1104 ml
<b>Deviation:</b>	0,9011 ml	0,7202 ml	0,7928 ml
<b>Precision (CV):</b>	0,4506 %	0,1200 %	0,0793 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

4

**Date:**

09.08.2024

**Balance:**

Sartorius -  
MSE14202S-  
OCE-D0

DE-M