

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 100 BL

Grad. cylinder, 1000:10 ml, blue grad., A, conf. cert.




Batch: 49-24

Serial:

Measuring Nr.: 71121

Batch size: 150

Volume:	200,000 ml	600,000 ml	1000,000 ml
Tolerance:	± 5,0000 ml	± 5,0000 ml	± 5,0000 ml
Average value:	200,3923 ml	599,7999 ml	999,6565 ml
Deviation:	0,7136 ml	0,8902 ml	0,7029 ml
Precision (CV):	0,3568 %	0,1484 %	0,0703 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	09.12.2024
Balance:	Sartorius - MSE14202S- OCE-D0

