

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 100

**Graduated cylinder, 1000:10 ml, A, conf. cert.**




**Batch:** 31-24

**Serial:**

**Measuring Nr.:** 68401

**Batch size:** 280

Volume:	200,000 ml	600,000 ml	1000,000 ml
<b>Tolerance:</b>	± 5,0000 ml	± 5,0000 ml	± 5,0000 ml
<b>Average value:</b>	198,6818 ml	599,1915 ml	999,9405 ml
<b>Deviation:</b>	0,7624 ml	0,6101 ml	0,5644 ml
<b>Precision (CV):</b>	0,3812 %	0,1017 %	0,0564 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	06.08.2024
<b>Balance:</b>	Sartorius - MSE14202S- OCE-D0

