

Certificate of performance



Article: 3 715 100

Graduated cylinder, 1000:10 ml, A, conf. cert.




Batch: 49-24

Serial:

Measuring Nr.: 71238

Batch size: 280

Volume:	200,000 ml	600,000 ml	1000,000 ml
Tolerance:	± 5,0000 ml	± 5,0000 ml	± 5,0000 ml
Average value:	200,0663 ml	599,6656 ml	999,5646 ml
Deviation:	0,6183 ml	0,7026 ml	1,0426 ml
Precision (CV):	0,3091 %	0,1171 %	0,1043 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	18.12.2024
Balance:	Sartorius - MSE14202S- OCE-D0

