

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 505 USP

**graduated pipettes 5 ml class AS USP KB**




**Batch:** 33-24

**Serial:**

**Measuring Nr.:** 68865

**Batch size:** 280

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
<b>Tolerance:</b>	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
<b>Average value:</b>	4,9926 ml	2,9944 ml	0,9975 ml
<b>Deviation:</b>	0,0056 ml	0,0038 ml	0,0034 ml
<b>Precision (CV):</b>	0,1114 %	0,1280 %	0,3386 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	26.08.2024
<b>Balance:</b>	Sartorius - LE244S-OCE

