

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 910 BL

**Grad. pipettes, 25:0,1 ml, blue grad., AS, conf. cert.**




**Batch:** 25-24

**Serial:**

**Measuring Nr.** 67909

**Batch size:** 279

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	24,9732 ml	14,9940 ml	5,0051 ml
Deviation:	0,0187 ml	0,0165 ml	0,0124 ml
Precision (CV):	0,0748 %	0,1098 %	0,2482 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	10.07.2024
Balance:	Sartorius - LE244S-OCE

